

# NMCP HIGH-LEVEL DISINFECTION LOG FOR TROPHON EPR Serial # \_\_\_\_\_

**LOCATION:** \_\_\_\_\_

## INSTRUCTIONS:

**SONEX CARTRIDGE REPLACEMENT:** Every 30 days and after each purging cycle. Document: **LOT #, EXP DATE, DATE LOADED & CHANGE DATE.**

**CHEMICAL INDICATOR :** A new chemical indicator is to be used with each cycle. Document **LOT BATCH # & EXP DATE.**

Date & time from machine	Cartridge lot #/exp. date	Cartridge load date	Cartridge change date	Chemical indicator lot #/exp. date	Trophon® EPR printed label record	Patient label
	L# EXP:			L# EXP:		
	L# EXP:			L# EXP:		
	L# EXP:			L# EXP:		
	L# EXP:			L# EXP:		
	L# EXP:			L# EXP:		

Effective: 24 Sep 2013

**\*\*\*\*\*If a cycle fails or the chemical indicator fails the color change test, reprocess the probe. If cycle fails again, remove Trophon® EPR from use and contact Trophon® representative for service. Document in next row “Trophon® EPR removed from service.” When Trophon® EPR returned to service, begin a new row. \*\*\*\*\***